


# Helping reduce your risk of IAD and CAUTI



## PrimaFit™ External Urine Management System for Females

### What is PrimaFit?

PrimaFit is designed to help manage urinary incontinence for females. It is an alternative solution to incontinence pads and diapers, which can lead to skin breakdown and incontinence-associated dermatitis (IAD). In some cases, your PrimaFit may even replace an indwelling catheter, which is the #1 risk factor for catheter-associated urinary tract infections (CAUTI).



### How does it protect my skin?

When properly positioned, PrimaFit absorbs urine and moves it away from your skin. The urine is suctioned into a collection canister, helping keep your skin dry.

Approximately  
**75%** of  
hospital-acquired UTIs

are associated with  
an indwelling  
urinary catheter.<sup>1</sup>

REFERENCES: 1. Centers for Disease Control and Prevention, Catheter-associated Urinary Tract Infections (CAUTIs). Available at: [https://www.cdc.gov/hai/ca\\_uti/uti.html](https://www.cdc.gov/hai/ca_uti/uti.html). Accessed August 24, 2017.

Sage, now part of Stryker

25615B © Sage Products LLC 2018

**stryker**


# Pour aider à réduire le risque de DI et d'IVU associée au cathéter



## Dispositif externe de prise en charge de l'urine PrimaFit™ pour femmes

### Qu'est-ce que PrimaFit?

PrimaFit est conçu pour simplifier les soins de l'incontinence urinaire chez la femme. Il offre une solution de rechange aux serviettes et culottes d'incontinence, qui peuvent causer une dégradation cutanée et une dermatite d'incontinence (DI). Dans certains cas, PrimaFit peut même remplacer la sonde à demeure, principal risque d'infection des voies urinaires (IVU) associée au cathéter.



### Comment protège-t-il la peau?

Correctement installé, PrimaFit absorbe l'urine et l'éloigne de la peau. L'urine est aspirée dans un récipient collecteur, pour mieux garder la peau au sec.

**RÉFÉRENCE** : 1. Centers for Disease Control and Prevention, Catheter-associated Urinary Tract Infections (CAUTIs). Disponible à [https://www.cdc.gov/hai/ca\\_uti/uti.html](https://www.cdc.gov/hai/ca_uti/uti.html). Consulté le 24 août 2017.

**Sage, maintenant partie de Stryker**

25615B © Sage Products LLC 2018

**Environ**  
**75 %** des IVU  
**nosocomiales**

sont associées au port d'une sonde urinaire à demeure<sup>1</sup>.